# LepiMag™ Blood Cell-free DNA 추출 키트

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# LepiMag™ Blood Cell-free DNA Extraction Kit

## Manual User Guide

Catalog Number: LPMG48

Revision: A.0

---

## Table of Contents

- Specimen preparation and storage method ..... 5
  - Sample for use ..... 5
  - Sample storage conditions ..... 5
    - Whole blood ..... 5
    - Plasma ..... 5
- Kit contents and storage methods..... 6
- Required materials not supplied..... 6
- Pre-experimental setup..... 7
- LepiMag™ Blood Cell-free DNA manual extraction preparation..... 7
- LepiMag™ Blood Cell-free DNA manual extraction methods ..... 8
  - Workflow ..... 8
  - Prepare cell-free plasma samples ..... 8
  - Lysis and Binding ..... 9
  - Washing with LepiMag™ Cell-Free DNA Wash Solution..... 10
  - Washing with 80% ethanol..... 10
  - Elution ..... 11
- Troubleshooting..... 12
- Safety Data Sheet ..... 13
  - Section1: Identification of the substance/mixture and of the company/undertaking ..... 13

- Section 2: Harzards identification ..... 14
- Section 3: Composition/information on ingredients..... 16
- Section 4: First aid measures..... 16
- Section 5: Firefighting measures ..... 17
- Section 6: Accidental release measures ..... 17
- Section 7: Handling and storage ..... 18
- **Section 8: Exposure controls and personal protection ..... 19**
- **Section 9: Physical and chemical properties ..... 21**
- **Section 10: Stability and reactivity ..... 22**
- **Section 11: Toxicology information ..... 22**
- **Section 12: Ecological information ..... 23**
- **Section 13: Disposal considerations ..... 24**
- **Section 14: Transport information ..... 24**
- **Section 15 Regulatory information..... 25**
- **Section 16: Other information ..... 29**

## Specimen preparation and storage method

### ■ Sample for use

- Plasma separated from human whole blood collected through venipuncture.

### ■ Sample storage conditions

#### ■ Whole blood:

- Collect blood in a blood collection tube containing cell-free DNA (cfDNA) preservative (ex: Roche tube, Cat. No. 07785674001 or 07832397001). Storage temperature and duration of collected blood should be followed according to manufacturer guidelines.
- Immediate plasma separation after blood collection is recommended. Whole blood can be stored at 4–8°C for up to 24 hours before plasma separation.

#### ■ Plasma:

- Perform cfDNA extraction within 2 hours after plasma separation.
- If conducting cfDNA extraction within 2 hours is not possible, store plasma at -70°C to -90°C for a maximum duration of 3 years.
- Frozen plasma can be thawed at room temperature (or at 4°C) and used immediately without re-freezing.

## Kit contents and storage methods

Reagent	Volume (mL)	Cat. No	Storage
LepiMag™ Cell-free DNA Magnetic Beads	1.5	A59467C	4°C <sup>1</sup>
LepiMag™ Cell-free DNA Proteinase K	1.5	A59469C	Room temperature
LepiMag™ Cell-free DNA SDS Solution	8	A59468C	Room temperature
LepiMag™ Cell-free DNA Elution Solution	5	A59466C	Room temperature
LepiMag™ Cell-free DNA Lysis/Binding Solution	125	A59463C	Room temperature
LepiMag™ Cell-free DNA Wash Solution	100	A59465C	Room temperature

<sup>1</sup> Shipped at room temperature but stored at 4 degrees upon receipt.

## Required materials not supplied

Equipment	
Item	Details
Micropipette	10 µL
	20 µL
	100 µL
	200 µL
	1000 µL
Low-temperature centrifuge	Blood collection tube use
Centrifuge	1.5 mL microtube tube use
	5 mL Eppendorf tube use
	15 mL conical tube use
	50 mL conical tube use
Vortexer	1.5 mL microtube use
Heating bath	
DynaMag™-2 Magnet	Invitrogen™ 12321D
DynaMag™-50 Magnet	Invitrogen™ 12302D

Consumables	
Item	Details
Conical tube	15 mL, 50 mL
Snap tube	5 mL
Sterile microtube	1.5 mL
Filtered tips	10 µL
	200 µL
	1000 µL
Cooling rack	1.5 mL microtube usage

Reagents	
Item	Details
Ethanol	80%

## Pre-experimental setup


- Inspect kit components and experimental equipment for contamination or damage before using the product.
- The operator is required to wear nitrile gloves, and the testing surface should be thoroughly cleaned using 70% ethanol before conducting experiments to eliminate foreign particles and contamination.

## LepiMag™ Blood Cell-free DNA manual extraction preparation

- Set the equipment to be used to the appropriate temperature according to the manufacturer’s manual.
- Perform all steps at room temperature (15~25°C) unless otherwise noted.
- When pipetting samples up and down, avoid creating bubbles.
- Vortex LepiMag™ Cell-free DNA Magnetic Beads to fully resuspend them before use.
- Incubate LepiMag™ Cell-free DNA Lysis/Binding Solution and LepiMag™ Cell-free DNA Wash Solution at 37°C for one hour if precipitates are visible. This can happen if storage temperatures are too low.

## LepiMag™ Blood Cell-free DNA manual extraction methods

### ■ Workflow

### ■ Prepare cell-free plasma samples

- Collect blood in a blood collection tube containing cfDNA preservative (Roche tube). Label barcode (patient number, age, sex, blood collection date, etc. can be directly written)
- Gently mix the blood collection tube on a roller mixer for 5 minutes or gently invert the tube 4-5 times.
  - ❖ **CAUTION:** Excessive inverting of blood samples can potentially cause DNA contamination during cfDNA extraction.
  - ❖ **CAUTION:** The collected whole blood should be handled in accordance with the storage temperature and duration recommended by the manufacturer of the blood collection tube.
- Centrifuge for 10 minutes at 1,600g - 2,500g, 2°C - 8°C using a low-temperature centrifuge.
- Transfer the plasma to a new 15 mL conical tube or 5 mL Eppendorf tube, being careful not to transfer the buffy coat formed in the middle of the blood collection tube. To prevent buffy coat contamination, it is recommended that only 90% of the plasma portion be transferred to the 15 mL conical tube or 5 mL Eppendorf tube.
- Centrifuge the 15 mL conical tube or 5 mL Eppendorf tube containing plasma for 10 minutes at 14,000g -

16,000g, 2°C - 8°C. If high-speed centrifugation is not possible, centrifuge for 20 minutes at 3,000g - 5,000g, 2°C - 8°C.

- Transfer only the supernatant to a new 5 mL snap tube. If immediate cfDNA extraction is not possible, store the plasma in a freezer at or below -70°C to -90°C for long term storage.

## ■ Lysis and Binding

- Add the following reagents to a 50 mL conical tube depending on the plasma volume according to Table 1.

[Table 1]

Reagents	Plasma	
	1.5 mL	3 mL
LepiMag™ Cell-free DNA Proteinase K	9 µL	18 µL
Plasma sample	1.5 mL	3 mL
LepiMag™ Cell-free DNA SDS Solution <sup>[1]</sup>	75 µL	150 µL

[1] Do not add SDS directly to the Proteinase K solution, to avoid inactivation of the Proteinase K.

- Mix thoroughly all reagents by vortexing for 5 seconds and incubate in the heating bath at 60°C for 20 minutes.
  - ❖ A thermomixer can be used as an alternative to the heating bath.
- After incubation is complete, place the tube containing the plasma sample on ice for 5 minutes to cool.
- Add the reagents according to Table 2 directly into 50 mL conical tube.

[Table 2]

Reagents	Plasma	
	1.5 mL	3 mL
LepiMag™ Cell-free DNA Lysis/Binding Solution	1.25 mL	2.5 mL
LepiMag™ Cell-free DNA Magnetic Beads	15 µL	30 µL

- ❖ **CAUTION:** Thoroughly vortex the LepiMag™ Cell-free DNA Magnetic Beads before use to avoid low extraction efficiency.
- Invert the tube 10 times, then vigorously vortex for 10 minutes.
  - ❖ **CAUTION:** Insufficient vortexing may result in low cfDNA yield

- Place the tube on a DynaMag™-50 Magnet (Invitrogen™ 12302D). To optimize the interaction between the magnetic beads and the magnet, tilt the magnetic rack/tube to increase the contact area between the solution and the magnet. Allow the setup to be still for at least 5 minutes, or until the magnetic beads and supernatant are completely separated.
- Once the magnetic beads and supernatant have completely separated, carefully discard the clear supernatant.

### ■ Washing with LepiMag™ Cell-Free DNA Wash Solution

- Using LepiMag™ Cell-free DNA Wash Solution, resuspend the magnetic beads according to the volume specified in Table 3.

[Table 3]

Reagents	Plasma	
	1.5 mL	3 mL
LepiMag™ Cell-free DNA Wash Solution	0.5 mL	1 mL

- Transfer the resuspended Magnetic Bead slurry to a new 1.5 mL microcentrifuge tube.
- Place the tube on the DynaMag™-2 Magnet. Wait for at least 2 minutes until the magnetic beads and supernatant are fully separated, then carefully remove the clear supernatant with a pipette.
- Remove the tube from the DynaMag™-2 Magnet and add LepiMag™ Cell Free DNA Wash Solution in the volume specified in Table 3. Vortex for 30 seconds and spin down.
- Place the tube back on the DynaMag™-2 Magnet. Wait for at least 2 minutes until the magnetic beads and supernatant are fully separated, carefully remove the clear supernatant with a pipette.

### ■ Washing with 80% ethanol

- Remove the tube from the DynaMag™-2 Magnet and add 1 mL of 80% ethanol along the tube wall. Vortex for 30 seconds and spin down.
- Place the tube on the DynaMag™-2 Magnet. Wait until the magnetic beads and supernatant are fully separated, carefully remove the clear supernatant with a pipette.

- Repeat the above 2 steps.
- Place the tube on the DynaMag™-2 Magnet, open the cap and dry for 5 minutes.
- Remove as much remaining liquid as possible from the tube.

## ■ Elution

- Add 50 µL of LepiMag™ Cell-free DNA Elution Solution to the sample in the 1.5 mL tube. Vortex for 5 minutes using a vortexer and spin down.
- Place the 1.5 mL tube on the DynaMag™-2 Magnet and wait for at least 2 minutes or until the magnetic beads and supernatant are fully separated.
- Carefully transfer the clear supernatant to a new, sterile 1.5 mL tube using a pipette.

## Troubleshooting

Observation	Possible cause	Recommended action
Lower yield than expected	The LepiMag™ Cell-free DNA Magnetic Beads were not properly stored	Remove the LepiMag™ Cell-free DNA Magnetic Beads from the kit and store them at 2–8°C. Do not freeze the beads. Allow the beads to warm to room temperature before use.
	An insufficient vortexing of LepiMag™ Cell-free DNA Magnetic Beads	Vortex the tube containing the magnetic beads thoroughly immediately before use.
		<p><b>Not thoroughly vortexed beads seen from bottom of tube</b>      <b>Thoroughly vortexed beads seen from bottom of tube</b></p> 
	The LepiMag™ Cell-free DNA Magnetic Beads are not optimally dried	Drying times may vary depending on the amount of beads used and the environment. Lower volumes of beads require less time to dry. Airflow and humidity in the immediate environment may shorten or lengthen the optimal bead drying time.
Overdried beads will stick to the wall of the microtube and be difficult to re-suspend.		
Underdried beads may carry ethanol into the eluate and negatively impact downstream applications.		
	The sample contains low amount of cfDNA	Increase the starting sample volume.
Magnetic bead carryover	Loose beads present in the eluate or inadvertently transferred	Be sure to leave the tube on the magnetic stand when removing the eluate containing the cfDNA.

## Safety Data Sheet

Safety Data Sheet according to Standards for Classification and Labeling of Chemical Substances and Material Safety Data Sheet (MoEL Public Notice No. 2023-9)

Submission Number: No data

---

### Section1: Identification of the substance/mixture and of the company/undertaking

#### **SECTION 1: Identification of the substance/mixture and of the company/undertaking**

##### Identification of the substance or mixture

Product code     A59463C

Product name     LepiMag™ Cell-free DNA Lysis/Binding Solution

##### Relevant identified uses of the substance or mixture and uses advised against

Relevant identified uses

For research use only.

Uses advised against

No information available.

##### Company/undertaking identification

Lepidyne Co., Ltd.

411 Seoul Forest A Tower, 26, Sangwon 1-gil,

Seongdong-gu, Seoul South Korea

Tel: + 82-2-2138-3818

Fax: + 82-2-2138-3828

## Section 2: Hazards identification

### SECTION 2: Hazards identification

#### GHS Classification

##### Health hazards

Acute oral toxicity	Category 4
Acute dermal toxicity	Category 1
Acute inhalation toxicity – gas	Category 1

##### Physical hazards

Not Hazardous

##### Environmental hazards

Chronic aquatic hazard	Category 3
------------------------	------------

#### Label elements

##### Signal Word

DANGER

##### Hazard pictograms

##### Hazard Statements

H302 – Harmful if swallowed

H314 – Causes severe skin burns and eye damage

H412 – Harmful to aquatic life with long lasting effects

#### Precautionary Statements

##### Prevention

P260 – Do not breathe dust/fume/gas/mist/vapours/spray

P264 – Wash hands thoroughly after handling

P270 – Do not eat, drink or smoke when using this product

P273 – Avoid release to the environment

P280 – Wear protective gloves/protective clothing/eye protection/face protection

##### Response

P303 + P361 + P353 – IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower

P301 + P310 – IF SWALLOWED: Immediately call a POISON CENTER or doctor/physician

P305 + P351 + P338 – IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing

P330 – Rinse mouth

P310 – Immediately call a POISON CENTER or doctor/physician

P304 + P340 – IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing

P301 + P330 + P331 – IF SWALLOWED: Rinse mouth. Do NOT induce vomiting

##### Storage

Not Applicable

##### Disposal

P501 – Dispose of contents/ container to an approved waste disposal plant

Other hazards

Contains a known or suspected endocrine disruptor.

### Section 3: Composition/information on ingredients

#### SECTION 3: Composition/information on ingredients

Chemical Name	CAS No	Common name	EINECS-No.	Weight-%
Thiocyanic acid, compound with guanidine(1:1)	593-84-0	Not applicable	209-812-1	35-45
Triton X-100	9002-93-1	Not applicable	-	16-24

Contact with acids or bleach liberates toxic gases. DO NOT ADD acids or bleach to any liquid wastes containing this product. We recommend handling all chemicals with caution.

### Section 4: First aid measures

#### SECTION 4: First aid measures

#### Description of first aid measures

Skin contact	Wash off immediately with plenty of water for at least 15 minutes. Remove and wash contaminated clothing and gloves, including the inside, before re-use. Immediate medical attention is required.
Eye contact	Rinse immediately with plenty of water, also under the eyelids, for at least 15 minutes. Immediate medical attention is required.
Ingestion	Never give anything by mouth to an unconscious person. Do not induce vomiting without medical advice. Get medical attention if symptoms occur.
Inhalation	Remove to fresh air. If not breathing, give artificial respiration. If symptoms persist, call a doctor.
Notes to Physician	Treat symptomatically.

Most important symptoms and effects, both acute and delayed

H302 – Harmful if swallowed

H314 – Causes severe skin burns and eye damage

H318 – Causes serious eye damage.

H412 – Harmful to aquatic life with long lasting effects

Indication of any immediate medical attention and special treatment needed

IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/ shower. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell.

## Section 5: Firefighting measures

### **SECTION 5: Firefighting measures**

Extinguishing media

Suitable extinguishing media

Water spray. Carbon dioxide (CO<sub>2</sub>). Foam. Dry chemical.

Unsuitable extinguishing media

No information available.

Special hazards arising from the substance or mixture

No data

Advice for firefighters

Wear self-contained breathing apparatus and protective suit.

## Section 6: Accidental release measures

### **SECTION 6: Accidental release measures**

Personal precautions, protective equipment and emergency procedures

Ensure adequate ventilation

Avoid contact with skin, eyes or clothing

Use personal protection equipment

See section 8 for more information

#### Environmental precautions

Prevent product from entering drains. Do not allow material to contaminate ground water system.

#### Methods and material for containment and cleaning up

Soak up with inert absorbent material.

#### Reference to other sections

See section 8 for more information.

### Section 7: Handling and storage

#### **SECTION 7: Handling and storage**

##### **Precautions for safe handling**

Always wear recommended Personal Protective Equipment. See section 8 for more information. Do not get in eyes, on skin, or on clothing. Do not ingest. If during normal use the material presents a respiratory hazard, use adequate ventilation and/or wear appropriate respirator.

##### **Conditions for safe storage, including any incompatibilities**

Keep in a dry, cool and well-ventilated place. Keep in properly labelled containers. Store in accordance with local regulations.

##### **Specific end use(s)**

For research use only.

#### **SECTION 8: Exposure controls and personal protection**

## Section 8: Exposure controls and personal protection

### Control parameters

Chemical Name	Australia STEL	Australia	Australia – Skin Absorbers	Australia – OES – Peak
Thiocyanic acid, compound with guanidine (1:1) 593-84-0	No data	No data	No data	No data
Triton X-100 9002-93-1	No data	No data	No data	No data

Chemical Name	New Zealand – Workplace Exposure Limits	New Zealand TWA	New Zealand STEL	New Zealand Ceiling Limits
Thiocyanic acid, compound with guanidine (1:1) 593-84-0	No data	No data	No data	No data
Triton X-100 9002-93-1	No data	No data	No data	No data

Chemical Name	China OEL (MAC)	China – Skin Absorbers	Hong Kong OEL
Thiocyanic acid, compound with guanidine (1:1) 593-84-0	No data	No data	No data
Triton X-100 9002-93-1	No data	No data	No data

Chemical Name	Korea OEL (CEILING)	Korea OEL (TWA)	Korea OEL (STEL)
Thiocyanic acid, compound with guanidine (1:1) 593-84-0	No data	No data	No data
Triton X-100 9002-93-1	No data	No data	No data

Chemical Name	Korea OEL (Skin Notations)	Korea PEL (STELs)	Korea PEL (TWAs)
Thiocyanic acid, compound with guanidine (1:1) 593-84-0	No data	No data	No data
Triton X-100 9002-93-1	No data	No data	No data

Chemical Name	Korea OEL (Carcinogens)	Korea OEL (Germ Cell Mutagenicity)	Korea OEL (Reproductive Toxicity)

Thiocyanic acid, compound with guanidine (1:1) 593-84-0	No data	No data	No data
Triton X-100 9002-93-1	No data	No data	No data

Chemical Name	Taiwan OEL – TWA	Taiwan OEL (Skin notations, Carcinogens)	Taiwan OEL (STELs)	Taiwan OEL (Ceilings)
Thiocyanic acid, compound with guanidine (1:1) 593-84-0	No data	No data	No data	No data
Triton X-100 9002-93-1	No data	No data	No data	No data

Chemical Name	Singapore OEL (PEL)	Singapore OEL (STEL)	Malaysia – OEL (TWA)	Malaysia – OEL (Ceilings)
Thiocyanic acid, compound with guanidine (1:1) 593-84-0	No data	No data	No data	No data
Triton X-100 9002-93-1	No data	No data	No data	No data

Chemical Name	ACGIH OEL (TWA)	ACGIH OEL (STEL)
Thiocyanic acid, compound with guanidine (1:1) 593-84-0	No data	No data
Triton X-100 9002-93-1	No data	No data

Engineering Measures                      Ensure adequate ventilation, especially in confined areas

**Exposure controls**

Personal Protective Equipment

Use PPE that has been certified by the Korea Occupational Safety and Health Agency.

Respiratory protection                      In case of insufficient ventilation wear respirators and components tested and approved under appropriate government standards.

Hand protection                              Wear suitable gloves Glove material: Compatible chemical-resistant gloves.

Eye protection	Tight sealing safety goggles
Skin and Body Protection	Wear laboratory coat for body protection
Hygiene Measures	Handle in accordance with good industrial hygiene and safety practice

## Section 9: Physical and chemical properties

### SECTION 9: Physical and chemical properties

#### Information on basic physical and chemical properties

Physical state	Liquid	
Colour	No data	
Odour	No data	
Molecular Weight	No data	
Melting point / melting range	°C No data	°F No data
Boiling point / boiling range	°C No data	°F No data
Flammability (solid, gas)	No data	
Lower explosion limit	No data	
Upper explosion limit	No data	
Flash point	°C No data	°F No data
Autoignition Temperature	°C No data	°F No data
Decomposition temperature	°C No data	°F No data
pH	No data	
Evaporation rate	No data	
Viscosity	No data	
Solubility	No data	
Partition coefficient	No data	
n-octanol/water	No data	
Vapour Pressure	No data	
Specific gravity	No data	
Vapour density	No data	

Explosive properties	No data
Oxidising properties	No data
Particle characteristics	No data

Other information

Information with regards to physical hazard classes

No information available

**Other safety characteristics**

No information available

**Section 10: Stability and reactivity**

**SECTION 10: Stability and reactivity**

Reactivity	None known.
Chemical stability	Stable under normal conditions.
Possibility of hazardous reactions	Hazardous reaction has not been reported.
Conditions to avoid	Contact with acids or bleach toxic gas. DO NOT ADD acids or bleach to any liquid wastes containing this product.
Incompatible materials	No dangerous reaction known under conditions of normal use.
Hazardous decomposition products	No known hazardous decomposition products.

**Section 11: Toxicology information**

**SECTION 11: Toxicological information**

**Information on toxicological effects**

Chemical Name	Oral LD50	Dermal LD50	Inhalation LC50
Thiocyanic acid, compound with guanidine (1:1)	593mg/kg	No data available	No data available
Triton X-100	= 1800 mg/kg (Rat)	No data available	No data available

Principal Routes of Exposure

<b>Acute toxicity</b>	Harmful if swallowed
<b>Skin corrosion/irritation</b>	Causes severe burns
<b>Serious eye damage/irritation</b>	Causes serious eye damage
<b>Respiratory or skin sensitization</b>	Data are conclusive but insufficient for classification
<b>(STOT) – single exposure</b>	Data are conclusive but insufficient for classification
<b>Specific target organ toxicity</b>	
<b>(STOT) – repeated exposure</b>	Data are conclusive but insufficient for classification
<b>Carcinogenicity</b>	
<b>Germ cell mutagenicity</b>	Data are conclusive but insufficient for classification
<b>Reproductive toxicity</b>	Data are conclusive but insufficient for classification
<b>Aspiration hazard</b>	Data are conclusive but insufficient for classification

## Section 12: Ecological information

### SECTION 12: Ecological information

#### Toxicity

Hazardous to the Aquatic Environment.

Chemical Name	Toxicity to algae	Toxicity to daphnia and other aquatic invertebrates	Toxicity to fish	Microtox Data	log Pow
Thiocyanic acid, compound with guanidine (1:1)	No data available	No data available	No data available	No data available	No data available
Triton X-100	No data available	EC50=26 mg/L (48 h)	LC50 = 8.9 mg/L 96H LC50 = 4.0 mg/l 96H (Pimephales promelus)	No data available	logPow-2.7

Mobility in soil	No information available
Persistence and degradability	No information available
Bioaccumulative potential	No information available
Results of PBT and vPvB assessment	No information available

Other adverse effects

Contains a known or suspected endocrine disruptor

### Section 13: Disposal considerations

#### SECTION 13: Disposal considerations

#### Waste treatment methods

The generation of waste should be avoided or minimized wherever possible. Empty containers or liners may retain some product residues. This material and its container must be disposed of in according to approved disposal technique. Disposal of this product, its solutions or of any by-products, shall comply with the requirements of all applicable local, regional or national/federal regulations. Do not empty into drains. Do not dispose of waste into sewer.

### Section 14: Transport information

#### SECTION 14: Transport information

#### IATA / ADR / DOT-US / IMDG

Classified as dangerous in the meaning of transport regulations

UN number or ID number	1760
UN proper shipping name	Corrosive liquid, n.o.s. (Guanidine Thiocyanate)
Transport hazard class(es)	8
Packing group	III

Environmental hazards

yes

Special precautions for user

Not Applicable

Maritime transport in bulk according to IMO instruments

Not Applicable.

## Section 15 Regulatory information

### SECTION 15: Regulatory information

#### International Inventories

Chemical Name	US TSCA	Canada (DSL)	NDSL	PICCS (Philippines)
Thiocyanic acid, compound with guanidine (1:1)	Listed	Listed	-	Listed
Triton X-100	Listed	Listed	-	Listed

Chemical Name	ENCS (Japan)	IECSC – China Inventory of Existing Chemical Substances	South Korea (KECL)	AICS (Australia)
Thiocyanic acid, compound with guanidine (1:1)	Listed	Listed	-	Listed
Triton X-100	Listed	Listed	Listed	Listed

Chemical Name	NZIoC – New Zealand Inventory of Chemicals	EINECS (European Union)	ELINCS (European List of Notified Chemical Substances)
Thiocyanic acid, compound with guanidine (1:1)	Listed	Listed	-
Triton X-100	Listed	Listed	-

#### UN (United Nations) – Chemical Weapons Convention (CWC)

Not Applicable

#### Korea

#### Occupational Safety and Health Act

Chemical Name	ISHA – Substances to be controlled – Metals	ISHA – Substances to be controlled – Organic Substances	ISHA – Substances to be controlled – Acids and bases	ISHA – Substances to be controlled – Gas-phase substances
Thiocyanic acid, compound with guanidine (1:1)	Not Applicable	Not Applicable	Not Applicable	Not Applicable
Triton X-100	Not Applicable	Not Applicable	Not Applicable	Not Applicable

Chemical Name	ISHA – Harmful Substances Prohibited for Manufacturing, Importing, Transferring, or Supplying	ISHA – Harmful Substances Requiring Permission	Special management materials
Thiocyanic acid, compound with guanidine(1:1)	Not Applicable	Not Applicable	Not Applicable
Triton X-100	Not Applicable	Not Applicable	Not Applicable

Chemical Name	Harmful agents subject to work environment monitoring – Organic compounds	Harmful agents subject to work environment monitoring – Metals	Harmful agents subject to work environment monitoring – Acids and alkalis	Harmful agents subject to work environment monitoring – Gas-phase substances	Harmful agents subject to work environment monitoring – Dusts
Thiocyanic acid, compound with guanidine (1:1)	Not Applicable	Not Applicable	Not Applicable	Not Applicable	Not Applicable
Triton X-100	Not Applicable	Not Applicable	Not Applicable	Not Applicable	Not Applicable

Chemical Name	Harmful agents subject to workers requiring health examination – Organic compounds	Harmful agents subject to workers requiring health examination – Metals	Harmful agents subject to workers requiring health examination – Acids and alkalis	Harmful agents subject to workers requiring health examination – Gas-phase substances	Harmful agents subject to workers requiring health examination – Dusts
Thiocyanic acid, compound with guanidine (1:1)	Not Applicable	Not Applicable	Not Applicable	Not Applicable	Not Applicable
Triton X-100	Not Applicable	Not Applicable	Not Applicable	Not Applicable	Not Applicable

Chemical Name	Korea OEL (CEILING)	Korea OEL (TWA)	Korea OEL (STEL)	Korea OEL (Skin Notations)
Thiocyanic acid, compound with guanidine (1:1)	Not Applicable	Not Applicable	Not Applicable	Not Applicable
Triton X-100	Not Applicable	Not Applicable	Not Applicable	Not Applicable

Chemical Name	Korea OEL (Carcinogens)	Korea OEL (Germ Cell Mutagenicity)	Korea OEL (Reproductive)	Korea PEL (STELs)	Korea PEL (TWAs)
---------------	-------------------------	------------------------------------	--------------------------	-------------------	------------------

			Toxicity)		
Thiocyanic acid, compound with guanidine (1:1)	Not Applicable	Not Applicable	Not Applicable	Not Applicable	Not Applicable
Triton X-100	Not Applicable	Not Applicable	Not Applicable	Not Applicable	Not Applicable

Chemical Name	ISHA – Harmful or Dangerous Substances Subject to Submission of Process Safety Reports – Storage	ISHA – Harmful or Dangerous Substances Subject to Submission of Process Safety Reports – Manufacturing or Handling
Thiocyanic acid, compound with guanidine (1:1)	Not Applicable	Not Applicable
Triton X-100	Not Applicable	Not Applicable

### Chemicals Control Act

Chemical Name	Toxic substance	Prohibited substance	Restricted substance
Thiocyanic acid, compound with guanidine(1:1)	Not Applicable	Not Applicable	Not Applicable
Triton X-100	Not Applicable	Not Applicable	Not Applicable

Chemical Name	Threshold Values (%)	Quantity limit for manufacture/use	Quantity limit for custody/storage
Thiocyanic acid, compound with guanidine (1:1)	Not Applicable	Not Applicable	Not Applicable
Triton X-100	Not Applicable	Not Applicable	Not Applicable

Chemical Name	Pollutant Release and Transfer Registers (PRTR) – Group I	Pollutant Release and Transfer Registers (PRTR) – Group II
Thiocyanic acid, compound with guanidine (1:1)	Not Applicable	Not Applicable
Triton X-100	Not Applicable	Not Applicable

### Act on Registration, Evaluation, etc. of Chemicals (K-REACH)

Chemical Name	Ministry of Environment – CMR risk	Ministry of Environment – Critically Controlled Substance	Existing substances subject to registration	Substances Exempted from Registration or Notification – Annex 1	Substances Exempted from Registration or Notification – Annex 2
Thiocyanic acid, compound with guanidine (1:1)	Not Applicable	Not Applicable	Not Applicable	Not Applicable	Not Applicable

Triton X-100	Not Applicable	EDCs	Not Applicable	Not Applicable	Not Applicable
--------------	----------------	------	----------------	----------------	----------------

Safety Control of Dangerous Substances Act

Chemical Name	Safety Control of Dangerous Substances Act
Thiocyanic acid, compound with guanidine (1:1)	Not Applicable
Triton X-100	Not Applicable

Wastes Control Act

Chemical Name	Wastes Control Act
Thiocyanic acid, compound with guanidine (1:1)	Not Applicable
Triton X-100	Not Applicable

Singapore

Chemical Name	Singapore – Poisons List
Thiocyanic acid, compound with guanidine (1:1)	Not Listed
Triton X-100	Not Listed

Taiwan

Chemical Name	Toxic Chemical Substances Control Act – Classification and Control Levels	Toxic Chemical Substances Control Act – Threshold Regulated Quantities	Hazard Prevention Standard for Specified Chemical Substances
Thiocyanic acid, compound with guanidine(1:1)	Not Listed	Not Listed	Not Listed
Triton X-100	Not Listed	Not Listed	Not Listed

China

Chemical Name	China – Catalogue of Occupational Hazard Factors	China – Dangerous Goods List	China – Catalog of Hypertoxic Chemicals	China – List of Dangerous Chemicals for Priority Management
Thiocyanic acid, compound with guanidine (1:1)	Not Listed	Not Listed	Not Listed	Not Listed
Triton X-100	Not Listed	Not Listed	Not Listed	Not Listed

## Section 16: Other information

### SECTION 16: Other information

Reason for revision	SDS sections updated.
Revision number	2
Revision date	23-Aug-2023
Date of first issue	29-Jun-2015

For research use only.

#### References

- ECHA: <http://echa.europa.eu/>
- TOXNET: <http://toxnet.nlm.nih.gov/>
- eChemPortal: <http://www.echemportal.org/>
- LOLI database: <https://www.chemadvisor.com/loli-database>

#### Abbreviations and acronyms

TWA – Time-Weighted Average

OELs – Occupational Exposure Limits

STEL – Short Term Exposure Limit

DSL/NDSL – Canadian Domestic Substances List/Non-Domestic Substances List

KECL – Korean Existing and Evaluated Chemical Substances

ENCS – Japan Existing and New Chemical Substances

IECSC – China Inventory of Existing Chemical Substances

PICCS – Philippines Inventory of Chemicals and Chemical Substances

AICS – Australian Inventory of Chemical Substances

NZIoC – New Zealand Inventory of Chemicals

EINECS/ELINCS – European Inventory of Existing Chemical Substances/European List of Notified Chemical Substances

TSCA – United States Toxic Substances Control Act Section 8(b) Inventory

CEPA – Canadian Environmental Protection Act

EPA – Environmental Protection Agency

OSHA – Occupational Safety and Health Administration of the US Department of Labour

IATA – International Air Transport Association

DOT – Department of Transportation

IMDG – International Maritime Dangerous Goods

ACGIH – American Conference of Governmental Industrial Hygienists

NIOSH – National Institute for Occupational Safety and Health

AIHA – American Industrial Hygiene Association

HMIS – Department of Defense Hazardous Materials Information System

NTP – National Toxicology Program

IARC – International Agency for Research on Cancer

"The above information was acquired by diligent search and/or investigation and the recommendations are based on prudent application of professional judgment. The information shall not be taken as being all inclusive and is to be used only as a guide. All materials and mixtures may present unknown hazards and should be used with caution.

Since the Company cannot control the actual methods, volumes, or conditions of use, the Company shall not be held liable for any damages or losses resulting from the handling or from contact with the product as described herein.

THE INFORMATION IN THIS SDS DOES NOT CONSTITUTE A WARRANTY, EXPRESSED OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURPOSE"

End of Safety Data Sheet